


Navigates Tortuous, Tightly Curved
Anatomy with
**Precision, Control, and
Refined Finesse**



iSail Edge Align™

Differentiation Strategy for Vessel Optimization

- **Engineered for Precision**
iSail Edge Align™ is a uniquely engineered solution for treatment with precision

- **Edge Align™ Technology**
iSail ensures controlled and targeted expansion of intent-to-treat-vessel segment in a precise manner

- **Precision in Complex Bifurcations**
With iSail Edge Align™, you get
 - **Precise control**
 - **Reduced risk of unintended vessel injury**
 - **Procedural efficiency**, especially in complex bifurcations



- **New Class of Balloon**
Optimised shoulder length - 0.6 mm
- **Mastering Complexity with the Seamless Power of 3T Technology**
Transition-Free Design
0.41 mm Tapered Profile for Lesion Entry
2 mm Tri-Layered Atraumatic Tip
- **Empowering Excellence in Cardiac Procedures**
 - FlexDura Balloon Technology**
Unique balloon formation with excellent wrap/re-wrap capabilities
 - NavPro Technology**
Smooth navigation and durable performance with enhanced pushability and kink resistance
- A **Precision Tool** for segmental and edge-level optimization, ideal in complex bifurcation or left main cases



iSail Edge Align™

Ordering Information

Balloon Diameter (mm)	Balloon Length (mm)/Item References				
	06	08	10	12	15
2.50	SNC-25006	SNC-25008	SNC-25010	SNC-25012	SNC-25015
2.75	SNC-27506	SNC-27508	SNC-27510	SNC-27512	SNC-27515
3.00	SNC-30006	SNC-30008	SNC-30010	SNC-30012	SNC-30015
3.50	SNC-35006	SNC-35008	SNC-35010	SNC-35012	SNC-35015

Device Technical Specifications

Features

iSail Edge Align™

Proximal Shaft	2.0F
Distal Shaft	2.7F
Catheter Working Length	140 cm
Tip length	2.00 mm (Ø2.50 – 3.00 mm); 2.50 mm (Ø3.50 mm)
Crossing Profile (Ø3.0 mm)	0.83 mm
Coating	Hydrophilic
Nominal Pressure (NP)	12 atm
Rated Burst Pressure (RBP)	22 atm
Balloon Diameters (mm)	2.50, 2.75, 3.00, 3.50
Balloon Lengths (mm)	8, 10, 12, 15
Guiding Catheter Compatibility	5F
Compatible Guidewire	0.014"



Head Office

Plot No 143-A 1, Bommasandra Industrial Area,
Hebbagodi Village, Attibele, Hobli, Anekal Taluk
Bommasandra Industrial Estate, Bengaluru, Bengaluru
South, Karnataka, India, 560099.
Telephone No.: +91-120-4193839
FAX: +91-120-4193839
Email: info@innvolution.com
Website: www.innnvolution.com

Manufacturing Unit

Plot No. 105, Solitaire Industrial Park, Village
Dehmikalan, The Sanganer, Jaipur,
Rajasthan, India 303 007.
Manufacturing License No: MFG/MD/2022/000680

International Office

Marina Bay Financial Centre Tower 1, Level 11, 8 Marina
Boulevard, Singapore 018 981.